

RH 2010

Hall measurement system

PhysTech GmbH

Hall, DLTS, Customized Physical Measurement Equipment

Egilbert Str.2, D-85368 Moosburg

Tel.: 08761 74633 / Fax: 08761 74634

e-mail: info@PhysTech.de / WWW.PhysTech.de


Theoretical aspects and nomenclature

All mathematics that are used in the software of the RH 2010 Hall measurement system, can be calculated straight forward from basic physics or is published in the 50th by van der Pauw (specific resistivity measurement) and Hall (Hall effect measurement).

All calculations and formulars are based on the following sample structur.

Sample structure.

Top view.



A square sample with 4 Ohmic contact at the edges.

The area of the contacts should be small compared to the total sample area (smaller that shown above). A circled sample with 4 contacts symetrically set on the boundary is also possible.

The specific resistivity ρ and the Hall coefficient R_H are evaluated from measured 4-point resistivities at different contact configurations and R_H at different magnetic fields (R_H).

The 4-point resistivities are always be measured using V/I curves with preset currents using a variable current source and a voltage measurement circuit with automatic amplification. The current source and the voltage measurement circuit can be set using a special relais switch box to every contacts of the sample.

The V/I curves are automatically measured using 2 to 40 V/I points (user definable). The 4-point resistivities are calculated from the slope of the V/I curves using linear regression. The regression section can be selected automatiacly or manually. It will be saved with the measured data.


In the Hall (and v.d. Pauw) software menus and printouts the following definitions are used.

- I (contactpair) : = current is set to defined contactpair
 e.g. I(1,2) : = Current is set to the contacts 1 and 2 of the sample.
- V (contactpair) : = Voltage is measured at the defined contact pair
 e.g. V(3,4) : = Voltage is measured at sample contacts 3 and 4


The resistivity is defined by: R (contactpair I / contactpair V)
 e.g. $R(1,2 / 3,4)$: current contactpair = 1,2, voltage contactpair = 3,4

If both contactpairs are the same (2-point measurement) only one is given for the resistivity definition. e.g. $R(1,2)$: 2-point resistivity between the contacts 1 and 2. Typical results of a 2-point and 4-point V/I measurement are shown below.

2-point measurement



4-point measurement



Using these definitions and the results of van der Pauw (Philips Res. Rep. 3 1958) the specific resistivity can be calculated out of two 4-point measurements:

$$\rho = C \cdot D \cdot F \cdot 0.5 \cdot (R(1,2/4,3) + R(2,3/1,4))$$

$$C = \pi / \ln(2)$$

F = Correction factor for unsymmetric samples with $R(1,2/4,3) \gg R(2,3/1,4)$ (see publication)

D = Thickness of the sample or the layer


The software displays also the sheet resistivity $\rho' = \rho$ at $D = 1\text{ cm}$

The correction factor F can not be given is a closed function. It is correctly calculated and put in a table as a function of the symmetry factor Sym:

$$\text{Sym} = R(1,2/4,3) / R(2,3/1,4)$$

All plots of van der Pauw results will show the values for the correction function and the symmetry factor.

Two different v.d.Pauw configurations are possible (one is used above). Both configurations are automatically measured and ρ is calculated by averaging ρ_1 and ρ_2 . Offset voltages (s. below, 1. plot) caused by contacts or thermoelectric effects do not effect the values for the v.d. Pauw resistances.



For Hall measurements a contact configuration is used, that gives current and voltage perpendicular to each other. R(1,3/2,4) or R(2,4/1,3) as measurement configurations are possible and can be selected for the measurement. For ideal samples (ideal square symmetry) the measured voltage at 0 magnetic field should be 0 independantly of the used current. This gives also 0 for the resistivity. Real samples will still give a voltage, the so called misalignement voltage. In some cases, if this misalignement voltage is much higher than the Hallvoltage, this voltage can cause big systematical errors. Therefore the RH 2010 Hallsystem contains a hardware compensatation of this misalignement voltage and gives much better results on asymeric samples as Hall systems without it. The compensation mode can be switched on by software. If the magnetic field is applied to the sample (perpendicular), the Hall voltage V_H can be measured at the voltage contacts. For the configuration above V_H is given by:

$$V_H = R_H * I * B * D^{-1}$$


- V_H = Hall voltage
- R_H = Hall coefficient
- I = Measurement current
- D = Thickness of the sample or layer
- B = Magnetic field

The Hall coefficient R_H contains the sample parameters (see below) and will be the result of a Hall measurement. Using the definitions above, we get:

$$R_H = D * (R(1,3/2,4), B = + B_{max}) - (R(1,3/2,4), B = - B_{max}) / (2 * B_{max})$$

The Hall measurements will be done always at magnetic fields symetrically to 0. B_{max} can be set by the user. The alternative contact configuration R(2,4/1,3, B) can also be selected by the user. Similar to the v.d.Pauw measurement the Hall resistivities $R(1,3/2,4,B)$ are measured by V/I curves and calculated by the slope, using linear regression.


A typical result of a Hall measurement is shown on the right.



The Hall coefficient R_H can also be measured using more than 2 B-fields by measuring the Hall resistance as a function of the magnetic field. R_H can be calculated then from the slope of the $R(1,3/2,4,B)$ versus B plot by:

$$R_H = D * \Delta R(1,3/2,4,B) / \Delta B$$

A typical result of this magnetic field dependence Hall measurement is shown below.



Each Hall resistance in the plot above is still measured by an V/I curve and calculated by automatic regression.

The measurement results Rho and R_H can now be used for calculation of the sample parameters as carrier concentration, carrier type and carrier mobility. Starting with:

$$Rho = (q * (\rho + n*b) * \mu_n)^{-1}$$

$$R_H = (\rho - n*b^2) / (q * (\rho + n*b)^2)$$

- q : electroncharge
- n : electronconcentration
- p : holeconcentration
- μ_n : electron mobility
- μ_p : hole mobility
- b : = μ_n / μ_p

For one carrier type material ($n \gg p$ or $p \gg n$) the following approximations are valid:

$$\text{Rho} = (q \cdot n \cdot \mu_n)^{-1} \quad \text{or} \quad \text{Rho} = (q \cdot p \cdot \mu_p)^{-1}$$

$$\text{R}_H = -(q \cdot n)^{-1} \quad \text{or} \quad \text{R}_H = +(q \cdot p)^{-1}$$


This leads to the following formulars, that are used for calculating the sample parameters out of the measurement results Rho and R_H .

$$n, p = (\text{R}_H \cdot q)^{-1} \quad - = n, + = p$$





$$\mu_n, \mu_p = \text{R}_H / \text{Rho} \quad - = \mu_n, + = \mu_p$$

The sign of R_H is a measurement result. The software automatically transfers it into n or p. These calculated values are also always displayed on a measurement result plot.


A complete combined v.d.Pauw and Hall measurement is shown below.
(A very symetrically sample with very good Ohmic contacts compared to that ones shown above)




A complete sample characterization including all 2-point, all 4-point, v.d. Pauw and Hall measurements are shown in the last two pages. The plots are as printed from our software. The first page shows a measurement @ two magnetic fields, the second one a measurement under varied magnetic field.





Name = @A_001.HTB
Comm = High Res amplifier
ID = LC_3
rcID = 1001
Date = 20.06.2000
Type = p-Si
Thick. = 250.00 um
Temp = 295.00 K
Field = 0.21 T
Rho = 3.94E+01 Ohm*cm
Rho' = 1.58E+03 Ohm
Sym = 5.01E+00
F = 8.21E-01
Rh = 1.16E+04 cm ³ /As
Rh' = 4.67E+05 cm ² /As
p = 5.34E+14 cm ⁻³
my = 2.96E+02 cm ² /Vs





Name = @A_001.HTB
 Comm = High Res amplifier
 ID = LC_3
 rcID = 1001
 Date = 20.06.2000
 Type = p-Si
 Thick. = 250.00 um
 Temp = 295.00 K
 Field = 0.00 T
 R(1,2) = 6.89E+03
 R(2,3) = 7.73E+03
 R(3,4) = 3.88E+03
 R(4,1) = 4.18E+03
 R(1,3) = 4.59E+03
 R(2,4) = 7.60E+03
 R(1,2/4,3) = 1.40E+02
 R(2,3/1,4) = 7.08E+02
 R(3,4/2,1) = 1.41E+02
 R(4,1/3,2) = 7.05E+02
 R(1,3/2,4) = 5.65E+02
 R(2,4/3,1) = -5.68E+02



Name = @A_001.HTB
 Comm = High Res amplifier
 ID = LC_3
 rcID = 1001
 Date = 20.06.2000
 Type = p-Si
 Thick. = 250.00 um
 Cont. = 1,3/2,4
 Temp = 295.00 K
 Rho = 3.94E+01 Ohm*cm
 Rho' = 1.58E+03 Ohm
 Sym = 5.01E+00
 F = 8.21E-01
 Rh = 1.16E+04 cm³/As
 Rh' = 4.64E+05 cm²/As
 p = 5.38E+14 cm⁻³
 my = 2.94E+02 cm²/Vs

